



# THE EFFECT OF MAGNETIC FIELD ON POSITRON RANGE AND SPATIAL RESOLUTION IN AN INTEGRATED WHOLE-BODY TIME-OF-FLIGHT PET/MRI SYSTEM

Shih-ying Huang, Dragana Savic, Jaewon Yang, Uttam Shrestha, Youngho Seo, *Senior Member IEEE*  
UCSF Physics Research Laboratory, Department of Radiology and Biomedical Imaging, University of California, San Francisco


## INTEGRATED PET/MRI SYSTEM



GE Investigational PET/MR




## A CHARGED PARTICLE MOVING THROUGH A MAGNETIC FIELD



Force parallel to the magnetic field:  $F_{\parallel} = 0$   
Force not parallel to the magnetic field:  $F = qv \sin \theta$

**QUESTION:** Does the magnetic field affect the PET/MRI system spatial resolution for long-range positron emitters?

## GEANT4 MONTE CARLO SIMULATIONS




GEANT4 Monte Carlo simulation version 4.9.6.p02	Recorded positron annihilation position
Applied uniform magnetic field in the z-dir of a water cylinder	Fitted the positron spread profile with bi-exponential function
Isotropic distribution of $^{18}\text{F}$ , $^{124}\text{I}$ , and $^{68}\text{Ga}$	Evaluate full-width at tenth maximum (FVTM) of the point spread function


## GEANT4 Simulated Positron Track of $^{68}\text{Ga}$ in Water



## $^{68}\text{Ga}$ Point Spread Profile with 3 T Magnetic Field in Water



## Positron range vs. Magnetic field




## What did the simulations tell us?

The magnetic field may reduce in-plane positron spread but slightly increase positron spread in the direction of the magnetic field


The effect of magnetic field on positron spread is more noticeable for long-range positron emitters such as  $^{68}\text{Ga}$ ,  $^{124}\text{I}$ , and  $^{82}\text{Rb}$

## EXPERIMENTAL IMAGING




	GE Investigational PET/MR (3 T)	GE PET/CT Discovery VCT
<b>PET Crystal Size</b>		
In-plane (mm)	3.95	4.7
Axial (mm)	5.3	6.3
<b>PET Detector Ring Diameter (cm)</b>	60	88
<b>3D PET Acquisition (min per bed/# of bed)</b>	3.0 / 1 bed	5.0 / 1 bed
<b>PET Reconstruction Method (GE research toolbox)</b>	3D Fourier-rebinned filtered backprojection	3D Fourier-rebinned filtered backprojection
<b>PET Reconstructed Image Voxel Size (mm<sup>3</sup>)</b>	1.37 x 1.37 x 2.78	1.37 x 1.37 x 3.27


## Image Analysis Workflow



## Line source in Air



## Line source in Water



## What did the phantom imaging show?

The ratio of  $\text{FWHM}_z$  to  $\text{FWHM}_{\text{Avg}(x,y)}$  suggests a small effect on system spatial resolution from the magnetic field for  $^{68}\text{Ga}$

The 3 T magnetic field has nonsignificant effect on PET/MRI system spatial resolution

## ACKNOWLEDGMENTS

- Tremendous help from Vahid Ravanfart (nuclear medicine technologist) for phantom imaging
- Helpful conversation with the GE research team

This work is supported by: 

- GE Healthcare
- National Cancer Institute R01 CA154561